


PENETREX®

**THE ONLY SELF-RESEALING
ORAL SYRINGE PENETRABLE
LIQUID SAMPLE CONTAINER**



Perfect for Use with Oral Syringes

Easily insert the oral syringe into the pre-punched Penetrex seal to aspirate a solution. The cap seals on the plastic syringe tip to ensure a clean, spill-free transfer every time!



Reseals After Multiple Transfers

The patented Penetrex container effectively reseals after multiple uses to **prevent leaks and contamination.**



Prevent Accidental Spills

The risk of accidental spills and loss of valuable solution is virtually eliminated with the use of Penetrex containers.



Penetrex for the Lab

Andwin offers Penetrex containers fully sealed and pierceable by standard pipette tips. The “always-closed” condition speeds procedures and **virtually eliminates the risk of cross-contamination** common with flip-top and screw-cap containers that are prone to splashing and spilling.



Put the benefits of Penetrex to work for your everyday operations. Penetrex features can help increase productivity, improve user safety and prevent cross-contamination.



Complete Penetrex Packaging Solutions

Andwin offers a full range of Penetrex child-resistant bottles and oral syringes for liquid medications. Custom design and production is also available.

Contact Andwin for more information:

(800) 497-3113



Andwin Scientific A Division of the Andwin Corporation

6636 Variel Ave. Woodland Hills, CA 91303 Tel. 818.999.2828 Fax: 818.999.0111 www.andwin.com
PENETREX® is a registered trademark of the Andwin Corp. ©2006 The Andwin Corp. Protected by patents and patents pending.



An ISO Certified Company

